

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

Multibacket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

PATIENT ROOM EQUIPMENT

PATIENT ROOM EQUIPMENT

JB Medico products for the Beds / Patient Room, support the staff's work functions, hygiene requirements, and safety.

Our suspension articles enable flexible installation in the patient room.

JB Medico equipment ensures and supports the highest possible patient safety.

Based on targeted attention to increased hygiene, in which we have undergone care and treatment procedures in cooperation with the staff.

Our products are from the outset intended to ensure that enough room is available for the staff to carry out treatment and care tasks properly.

By designing our equipment so that they can, adapted to the needs and wishes of each department, one can meet one of the staff's most important professional requirements for the patient room, namely the ability to offer individual care and treatment, with the patient at the center

PATIENT ROOM EQUIPMENT

See the rest of the selection in the menu on the left side

CONTACT US FOR ADVICE AND OFFER ON HOSPITAL EQUIPMENT

If you would like advice or want to hear more about our hospital equipment, please call +45 7020 5826 or send us an email: info@jbmedico.com

Medical rail clamps

Medical rail clamps

For attachment of hospital equipment, medical rail clamps, used for the equipment rail located on the wall or bedside panel.

On the rail, different types of hospital equipment are attached using clamps, which provides fast, secure and flexible attachment.

There are several different standards of claws, and the different variants are often incompatible with each other.

When designing and ordering it is important to know which system to use.

JB-Medico makes available several different variants

JB-Medico has also developed claws 10 × 25 mm and 10 × 30 mm, which can handle most types of globes available on the market.

Can be mounted without tools and can be washed and disinfected in the dishwasher, as required

Together with leading specialists, we ensure our customers the most qualitative, hygienic and flexible product programs for working environments in hospitals and clinics, etc.

Made of necessity - designed to last.

We hope you can enjoy our products

CONTACT US FOR ADVICE AND OFFER ON HOSPITAL EQUIPMENT

Call phone +45-7020 5826 or send us an email at info@jbmedico.com

Multibracket, 2xØ6 holes, fixing device Ø18mm hole, JB 158-03-18

MOQ 2 PCS.

Download: [JB 158-03-18](#)

Multibracket with 2xØ6 clearance holes / Fit from 16-41mm. With Fixing device, Ø18mm hole.

Outer specifications:

H: 36mm W: 51mm D: 134/96mm / Weight 460 g

Tree bush sizes can be purchased and handles can be purchased

1. JB 128-00-09: Plastic bush with hole size Ø9.6mm
2. JB 128-00-01: Plastic bush with hole size Ø13mm
3. JB 127-00-01: Plastic bush with hole size Ø16mm
4. JB 94-00-00: Adjustable Handle M6x20mm
5. JB 163-00-01: A T-slot bracket for multi-bracket

Made of necessity - designed to last.

CONTACT US FOR ADVICE AND OFFER ON HOSPITAL EQUIPMENT

Call phone +45 7020 5826 or send us an email at info@jbmedico.com