

4-Position Transducer Mount, JB 289-00-00

	Description	P/N
	4-Position Transducer Mount Easy to clean and fits most transducers. Specifications: H: 90 mm / L: 185 mm / W: 0,4 kg Material: Aluminum Accessories: Multibracket with a T-slot bracket. Slide clamp half, locked with 1 ball clasp & T-slot DK Rail clamp, half with 1 ball clasp and T-slot EU T-Slot Wall bracket	JB 289-00-00 JB 158-00-163 JB 120-00-00 JB 146-00-00 JB 47-00-00

