

IV Hook stainless steel for infusion bags and IV fluids

Description	P/N
Drop bar for infusion bags, IV fluids with end cap. The rod fits all our existing fittings with the control for T-tracks Room for 5-6 drip bags / Easy to clean Specifications: Tube: Ø8 mm L: 200 mm W: 0,1 kg Material: Electro-polished Stainless steel Accessories: Multibracket with a T-slot bracket. JB 158-00-163 or Slide clamp half, with one ball clasp and T-slot JB 146-00-00	JB 205-08-05

