

IV- Stands/IV-Poles

P/N JB 320-00-02	P/N JB 321-00-02	P/N JB 317-00-05	P/N JB 319-00-03	P/N JB 317-00-13
Spec.	Spec.	Spec.	Spec.	Spec.
Bended, long outer tube L: 1281mm D: Ø32mm SST AISI 304	Bended, short outer tube L: 1076mm D: Ø32mm SST AISI 304	Straight outer tube, long L: 1281mm D: Ø32mm SST AISI 304	Straight outer tube, short L: 1076mm D: Ø32mm SST AISI 304	Telescopic pole fits IV – Stands L: 940mm D: Ø20mm SST AISI 304
				

