

STAINLESS STEEL ELECTRIC ADJUSTABLE TABLES

Instrument and Mayo tables

	Description	P/N
	Arm table, JB 100 30 The Remeda arm table is specially developed and suited for shoulder operations. Height adjustment: From 720mm to 1220mm (Electrical) Table top size: 248mm x 600mm (LxW)	JB 100 30
	Mayo table 3 legs, JB 100 21 Assistance Table is specifically designed to simplify the work in the OR room. Variants Height/JUST Table top size Lift 900-1400mm. 500x400mm (30kg) 900-1400mm 600x400mm (30kg) 900-1400mm 700x500mm (50kg) Swivel table top 900-1400mm 600x400mm (30kg) 900-1400mm 700x500mm (30kg)	JB 100 21-31 JB 100 21-32 JB 100 21-33 JB 100 21-42 JB 100 21-43
	Instrument table, two legs JB 100 22 The table tops are more significant than the assistance table and can carry more weight thanks to the solid frame Variants Height/JUST Table top size Lift 800-1200mm. 800x520mm (50kg) 900-1200mm 900x550mm (50kg) 900-1200mm 1070x570mm (50kg) 900-1400mm 800x520mm (50kg) 900-1400mm 900x550mm (50kg) 900-1400mm 1070x570mm (50kg)	JB 100 22 31 JB 100 22-32 JB 100 22-33 JB 100 22-41 JB 100 22-42 JB 100 22-43

Instrument and Mayo tables continued

	Description	P/N
 	Instrument table, JB 100 27 The 2 tier theatre instrument/back table is specially developed to simplify work in the operating Variants Height/JUST Table top size Lift 720-970mm. 1200x750mm (160kg) Placing shelf 1200x550mm 720-970mm. 1500x750mm (160kg) Draws shelf 1500x550mm 720-970mm. 1800x750mm (160kg) Draws shelf 1800x550mm Lower shelf 800x520mm (30kg) 900x550mm (50kg) 1070x570mm (30kg)	JB 100 27-01 JB 100 27-02 JB 100 27-03 JB 100 24-U3 JB 100 24-U4 JB 100 24-U6
	Instrument tables, two legs, JB 100 23 Assistance table is specially designed to simplify the work in the operating theatre. This table has the same features as the JB 100 20 table except that the pillar moved to the back of the tabletop. Variants Height/JUST Table top size Lift 800-1200mm. 800x520mm (50kg) 800-1200mm. 900x550mm (50kg) 800-1200mm. 1070x570mm (50kg) 800-1200mm. 1200x400mm (50kg) 900-1400mm. 800x520mm (50kg) 900-1400mm. 900x550mm (50kg) 900-1400mm. 1070x570mm (50kg) 900-1400mm. 1200x570mm (50kg)	JB 100-23-31 JB 100-23-32 JB 100-23-33 JB 100-23-34 JB 100-23-41 JB 100-23-42 JB 100-23-43 JB 100-23-44

Instrument and Mayo tables, continued

	Description	P/N
	Instrument tables, four legs, JB 100 20 The height can easily be adjusted using two buttons under the tabletop. Variants Height/JUST Table top size Lift 800-1200mm 800x520mm (50kg) 800-1200mm 900x550mm (50kg) 800-1200mm 1070x570mm (50kg) 800-1200mm 1200x570mm (50kg) 800-1400mm 800x520mm (50kg) 800-1400mm 900x550mm (50kg) 800-1400mm 1070x570mm (50kg) 800-1400mm 1200x570mm (50kg)	JB 100 20-31 JB 100 20-32 JB 100 20-33 JB 100 20-34 JB 100 20-41 JB 100 20-42 JB 100 20-43 JB 100 20-44
	Instrument table, heavy loads, 80 kg, two legs, JB 100 24 Has capacity for higher weight than standard instrument tables, and is thus particularly suitable for major operational treatments. For example, before orthopaedics. Surgery Variants Height/JUST Table top size Lift 800-1200mm 1070x570mm (80kg) 800-1200mm 1200x570mm (80kg) 800-1200mm 1200x900mm (80kg) 800-1200mm 1400x600mm (80kg) 900-1400mm. 1070x570mm (80kg) 900-1400mm. 1200x570mm (80kg) 900-1400mm. 1200x900mm (80kg) 900-1400mm. 1400x600mm (80kg) Lower shelf Table top size 800x520mm (30kg) 900x550mm (30kg) 1070x570mm (30kg)	JB 100-24-33 JB 100-24-34 JB 100-24-36 JB 100-24-38 JB 100-24-43 JB 100-24-44 JB 100-24-46 JB 100-24-48 JB 100 24-U3 JB 100 24-U4 JB 100 24-U6

Instrument and Mayo tables continued

Description: Instrument table with folded side tables JB 100 24-XX-FT			P/N
Instrument board with folded side tables is specially designed to simplify work in the operating theatre.			
The height can easily be adjusted using two buttons under the tabletop.			
Variants			
Height	JUST Table top size	Lift	
800-1200mm	1070 x 570mm + 2 x 400x570mm (15x80x15kg)	JB 100 24-33-FT	
800-1200mm	1200 x 570mm + 2 x 400x570mm (15x80x15kg)	JB 100 24-34-FT	
800-1200mm	1400 x 600mm + 2 x 400x600mm (15x80x15kg)	JB 100 24-38-FT	
900-1400mm	1070 x 570mm + 2 x 400x570mm (15x80x15kg)	JB 100 24-43-FT	
900-1400mm	1200 x 570mm + 2 x 400x570mm (15x80x15kg)	JB 100 24-44-FT	
900-1400mm	1400 x 600mm + 2 x 400x600mm (15x80x15kg)	JB 100 24-48-FT	



