






Multibracket, turned upside down, fit from 16-41mm

	Description Multibracket can be turned upside down, fit from 16-41mm Outer specifications: H: 36mm W: 51mm D: 134/96mm Weight: 230 g	P/N JB 158-00-00 MOQ 4 PCS.
	Multibracket with 2 x Ø6 clearance holes That allows you to put together div. Equipment, DIN rail and columns Outer specifications: H: 36mm W: 51mm D: 134/96mm Weight: 230 g	JB 158-00-05 MOQ 4 PCS.
	Double Multibracket for mounting, e.g. bed IV pole There are two pieces included. Outer specifications: Weight: 460 g Material: Aluminum, RF steel & plastic	JB 158-06-05 MOQ 2 PCS.
	Multibracket with a T-slot bracket can be turned upside down, fit from 16-41mm. Outer specifications: H: 36mm W: 51mm D: 134/96mm Weight: 230 g	JB 158-00-163 MOQ 4 PCS.
	Multibracket with 2xØ6 clearance holes Fit from 16-41mm. With Fixing device, Ø18mm hole. 3 bush sizes can be purchased Handles can be purchased Outer specifications: H: 36mm W: 51mm D: 134/96mm Weight 460 g	JB 158-03-18 MOQ 2 PCS.

Multibracket, turned upside down, fit from 16-41mm


Universal clamp for medical rails, round and square tubes, with eight drill holes on each side.

JB Medico develops the multibracket in designing a universal bracket suitable for mounting on various medical rails, round and square tubes.

The bracket is in anodised aluminum with JB-threaded holes, which fits a variety of JB-Medico's other products.

Multibracket screwed on a rail. Clamped on columns, tables for hypodermic needles, on DIN EU 10x25mm rail profiles and glovebox holders etc.

Can be mounted both in horizontal and vertical position.



FIXING DEVICE FOR SLIDE CLAMPS AND MULTI BRACKETS, Ø25MM

P/N: JB 205-00-00

Aluminum (anodised)

Prepared for suspension

Weight: 0,4 kg

Ø25mm Through hole without brass bushing

Specifications:

H: 61.5mm W: 45mm D: 35mm

Without brass bushing and adjustable stainless handle


Must be purchased as add-on products

Brass bushing: Ø25 x 20mm

P/N: JB 22-00-05

Stainless adjustable handle: M6 x 20mm

P/N: JB 94-00-00



Fixing device for slide clamps and multi-brackets, Ø18mm

P/N: JB 205-00-18





Aluminum (anodized)


Weight: 0,4 kg

Ø18mm hole



Accessories

	Description	P/N
	Plastic bush with hole size Ø9,6mm	JB 128-00-09
	Plastic bush with hole size Ø13mm	JB 128-00-01
	Plastic bush with hole size Ø16mm	JB 127-00-01
	Stainless adjustable handle M6 x 20mm	JB 94-00-00



GALLERY FIXING DEVICE FOR SLIDE CLAMPS AND MULTI BRACKETS Ø25 & Ø18mm



